

TdT Rabbit mAb


ABZ3192

Product Overview


Name	TdT Rabbit mAb
Accession(Primary)	P04053
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 59 kDa; Observed MW: 59 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:100-1:200; FC-1:20; IP-1:20

Target information (P04053)

Gene Name	DNTT
About Gene Symbol	



Western blot analysis of TdT in MOLT4 lysates using DNA Nucleotidyltransferase antibody.



Western blot analysis of TdT in 3T3 lysates using TdT antibody.