


Tetranectin Rabbit mAb

ABZ3193


Product Overview

Name	Tetranectin Rabbit mAb
Accession(Primary)	P05452
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 23 kDa; Observed MW: 23 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000


Target information (P05452)

Gene Name CLEC3B**About Gene Symbol**

Western blot analysis of Tetranectin in mouse lung lysates using Tetranectin antibody.



Western blot analysis of Tetranectin in MCF-7 lysates using Tetranectin antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Tetranectin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.