

TFE3 Rabbit mAb


ABZ3194

Product Overview

Name	TFE3 Rabbit mAb
Accession(Primary)	P19532
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 62 kDa; Observed MW: 62 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:100-1:200

Target information (P19532)

Gene Name	TFE3
About Gene Symbol	TFE3 (Transcription Factor Binding To IGHM Enhancer 3) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of TFE3 in HeLa, 293T, Raw264.7, Jurkat lysates using Transcription Factor E3 antibody.