

Thrombomodulin Rabbit mAb


ABZ3195

Product Overview


Name	Thrombomodulin Rabbit mAb
Accession(Primary)	P15306
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse, Rat
Calculated MW	Calculated MW: 62 kDa; Observed MW: 95 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:1000; IP-1:50-1:100

Target information (P15306)

Gene Name	Thbd
About Gene Symbol	



Western blot analysis of Thrombomodulin in rat lung lysates using Thrombomodulin antibody.



Western blot analysis of Thrombomodulin in 3T3 lysates using Thrombomodulin antibody.