

TIMELESS Rabbit mAb


ABZ3198

Product Overview

Name	TIMELESS Rabbit mAb
Accession(Primary)	Q9UNS1
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 139 kDa; Observed MW: 150 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:500; ICC/IF-1:20-1:50

Target information (Q9UNS1)

Gene Name	TIMELESS
About Gene Symbol	TIMELESS (Timeless Circadian Regulator) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Timeless in MCF-7, HeLa lysates using TIMELESS antibody.