

## TIMM44 Rabbit mAb


ABZ3199

### Product Overview


|                             |                                                                       |
|-----------------------------|-----------------------------------------------------------------------|
| <b>Name</b>                 | TIMM44 Rabbit mAb                                                     |
| <b>Accession(Primary)</b>   | O43615                                                                |
| <b>Host Species</b>         | Rabbit                                                                |
| <b>Clonality</b>            | Monoclonal                                                            |
| <b>Isotype</b>              | IgG                                                                   |
| <b>Species Reactivity</b>   | Human, Rat                                                            |
| <b>Calculated MW</b>        | Calculated MW: 51 kDa; Observed MW: 45 kDa                            |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; ICC/IF-1:200-1:500;<br>FC-1:200-1:500; IP-1:50-1:100 |

### Target information (O43615)

|                          |        |
|--------------------------|--------|
| <b>Gene Name</b>         | TIMM44 |
| <b>About Gene Symbol</b> |        |



Western blot analysis of TIMM44 in HeLa lysates using TIMM44 antibody.



Immunocytochemistry analysis of TIMM44(green) in HeLa using TIMM44 antibody, and DAPI(blue)