

TMEM192 Rabbit mAb

ABZ3200


Product Overview

Name	TMEM192 Rabbit mAb
Accession(Primary)	Q8IY95
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 31 kDa; Observed MW: 31 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:100-1:200

Target information (Q8IY95)

Gene Name TMEM192

About Gene Symbol



Western blot analysis of TMEM192 in K562 lysates using Transmembrane Protein 192 antibody.

