

Transmembrane Protein 30B Rabbit mAb

ABZ3201


Product Overview

Name	Transmembrane Protein 30B Rabbit mAb
Accession(Primary)	Q3MIR4
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 39 kDa; Observed MW: 39 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q3MIR4)

Gene Name TMEM30B

About Gene Symbol



Western blot analysis of TMEM30B in MCF-7 lysates using Transmembrane Protein 30B antibody.

