

TNFAIP3 Rabbit mAb

ABZ3202

Product Overview


Name	TNFAIP3 Rabbit mAb
Accession(Primary)	P21580
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 90 kDa; Observed MW: 82 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:50-1:100

Target information (P21580)


Gene Name	TNFAIP3
About Gene Symbol	



Immunocytochemistry analysis of TNFAIP3 (green) in HL-60 using TNFAIP3 antibody, and DAPI(blue).



Western blot analysis of TNFAIP3 in HeLa lysates using TNFAIP3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer tissue using TNFAIP3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.