


Tensin 4 Rabbit mAb

ABZ3203


Product Overview

Name	Tensin 4 Rabbit mAb
Accession(Primary)	Q8IZW8
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 77 kDa; Observed MW: 90 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q8IZW8)

Gene Name TNS4**About Gene Symbol**

Western blot analysis of TNS4/CTEN in Raw264.7, C6 lysates using Tensin 4 antibody.



Western blot analysis of TNS4/CTEN in K562
lysates using TNS4/CTEN antibody.