


TPP1 Rabbit mAb

ABZ3205

Product Overview

Name	TPP1 Rabbit mAb
Accession(Primary)	O14773
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 61 kDa; Observed MW: 61,48 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100

Target information (O14773)

Gene Name TPP1**About Gene Symbol**

Western blot analysis of TPP1 in HeLa, Jurkat, HepG2 lysates using TPP1 antibody.

