

Transglutaminase 2 Rabbit mAb


ABZ3206

Product Overview


Name	Transglutaminase 2 Rabbit mAb
Accession(Primary)	P21980
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 77 kDa; Observed MW: 77 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:2000-1:5000

Target information (P21980)


Gene Name	TGM2
About Gene Symbol	TGM2 (Transglutaminase 2) is a protease enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic



Immunohistochemistry analysis of paraffin-embedded mouse colon using Transglutaminase 2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Transglutaminase 2 in HUVEC lysates using Transglutaminase 2 antibody.



Western blot analysis of Transglutaminase 2 in rat
Brain lysates using Transglutaminase 2 antibody