

Translin Associated Protein X Rabbit mAb


ABZ3207

Product Overview


Name	Translin Associated Protein X Rabbit mAb
Accession(Primary)	Q99598
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 33 kDa; Observed MW: 33 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:50-1:100; IP-1:20-1:50

Target information (Q99598)

Gene Name	TSNAX
About Gene Symbol	TSNAX (Translin Associated Factor X) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of TRAX in 3T3, mouse testis lysates using Translin Associated Protein X antibody.



Western blot analysis of TRAX in K562, rat Brain, C6, HeLa lysates using TRAX antibody.