

## TARBP2 Rabbit mAb


ABZ3208

### Product Overview


|                             |                                                                                       |
|-----------------------------|---------------------------------------------------------------------------------------|
| <b>Name</b>                 | TARBP2 Rabbit mAb                                                                     |
| <b>Accession(Primary)</b>   | Q15633                                                                                |
| <b>Host Species</b>         | Rabbit                                                                                |
| <b>Clonality</b>            | Monoclonal                                                                            |
| <b>Isotype</b>              | IgG                                                                                   |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                                     |
| <b>Calculated MW</b>        | Calculated MW: 39 kDa; Observed MW: 39 kDa                                            |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:100-1:200;<br>ICC/IF-1:20-1:50; FC-1:50-1:100; IP-1:20-1:50 |

### Target information (Q15633)

|                          |                                                                                                                                                                                                                                                                 |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | TARBP2                                                                                                                                                                                                                                                          |
| <b>About Gene Symbol</b> | TARBP2 (TARBP2 Subunit Of RISC Loading Complex) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in infectious disease and is widely studied to elucidate cellular mechanisms and i |



Western blot analysis of TRBP in mouse liver lysates using TARBP2 antibody.



Western blot analysis of TRBP in HeLa, C6, MCF-7 lysates using TRBP antibody.