


TRIM21 Rabbit mAb

ABZ3209

Product Overview

Name	TRIM21 Rabbit mAb
Accession(Primary)	P19474
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 54 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100

Target information (P19474)

Gene Name TRIM21**About Gene Symbol**

Western blot analysis of TRIM21 in K562, HeLa lysates using TRIM21 antibody.

