

TRIM23/ARD1 Rabbit mAb

ABZ3210


Product Overview

Name	TRIM23/ARD1 Rabbit mAb
Accession(Primary)	P36406
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 64 kDa; Observed MW: 64 kDa
Recommended Dilution	WB-1:2000-1:5000

Target information (P36406)

Gene Name TRIM23

About Gene Symbol



Western blot analysis of TRIM23 in K562, C6, 3T3, HeLa lysates using TRIM23 antibody.

