

TXNL1 Rabbit mAb

ABZ3215


Product Overview

Name	TXNL1 Rabbit mAb
Accession(Primary)	O43396
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 32 kDa; Observed MW: 32 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200


Target information (O43396)

Gene Name TXNL1


About Gene Symbol



Western blot analysis of TXNL1 in K562, rat Brain, C6, 3T3, HeLa lysates using Thio redoxin Like Protein 1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using TXNL1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of TXNL1 (green) in HeLa using TXNL1 antibody, and DAPI (blue)

