

Tyrosine Hydroxylase Rabbit mAb


ABZ3216

Product Overview

Name	Tyrosine Hydroxylase Rabbit mAb
Accession(Primary)	P07101
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 59 kDa; Observed MW: 59 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:100-1:200

Target information (P07101)

Gene Name	TH
About Gene Symbol	



Western blot analysis of Tyrosine Hydroxylase in PC-12 lysates using Tyrosine 3 Monooxygenase antibody.