

U1C Rabbit mAb


ABZ3217

Product Overview


Name	U1C Rabbit mAb
Accession(Primary)	P09234
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 17 kDa; Observed MW: 21 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:200-1:500; IP-1:20-1:50

Target information (P09234)

Gene Name	SNRPC
About Gene Symbol	SNRPC (Small Nuclear Ribonucleoprotein Polypeptide C) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of U1C in rat lung lysates using U1C antibody.



Western blot analysis of U1C in THP-1, Jurkat, Hela, K562 lysates using U1C antibody.