

Ubiquitin Rabbit mAb


ABZ3218

Product Overview


Name	Ubiquitin Rabbit mAb
Accession(Primary)	P0CG47
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 26 kDa; Observed MW: 8 kDa
Recommended Dilution	WB-1:1000-1:20000; IHC-P-1:100-1:200; ICC/IF-1:200-1:500; FC-1:100-1:1000

Target information (P0CG47)


Gene Name	UBB
About Gene Symbol	



Immunocytochemistry analysis of Ubiquitin (green) in CEM using Ubiquitin antibody, and DAPI(blue).



Western blot analysis of Ubiquitin in rat brain, rat kidney lysates using Ubiquitin antibody.



Western blot analysis of Ubiquitin in 293 lysates
using Ubiquitin antibody