

## Ubiquitin Rabbit mAb

ABZ3219


### Product Overview

|                             |                                           |
|-----------------------------|-------------------------------------------|
| <b>Name</b>                 | Ubiquitin Rabbit mAb                      |
| <b>Accession(Primary)</b>   | P0CG47                                    |
| <b>Host Species</b>         | Rabbit                                    |
| <b>Clonality</b>            | Monoclonal                                |
| <b>Isotype</b>              | IgG                                       |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                         |
| <b>Calculated MW</b>        | Calculated MW: 26 kDa; Observed MW: 8 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:200-1:500       |


### Target information (P0CG47)

**Gene Name** UBB

**About Gene Symbol**



Western blot analysis of Ubiquitin in 293 lysates using Ubiquitin antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Ubiquitin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.