

## HAUS7 Rabbit mAb


ABZ3220

### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | HAUS7 Rabbit mAb                           |
| <b>Accession(Primary)</b>   | Q99871                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human                                      |
| <b>Calculated MW</b>        | Calculated MW: 40 kDa; Observed MW: 40 kDa |
| <b>Recommended Dilution</b> | WB-1:2000-1:5000                           |

### Target information (Q99871)

|                          |                                                                                                                                                                                                                                                                                                    |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | HAUS7                                                                                                                                                                                                                                                                                              |
| <b>About Gene Symbol</b> | HAUS7 (HAUS Augmin Like Complex Subunit 7) is a metabolic enzyme involved in ubiquitin-proteasome pathway. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of UCHL5IP in K562 lysates using HAUS7 antibody.