

Complex III Subunit 7 Rabbit mAb


ABZ3222

Product Overview


Name	Complex III Subunit 7 Rabbit mAb
Accession(Primary)	P14927
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 14 kDa; Observed MW: 14 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:100-1:200

Target information (P14927)

Gene Name	UQCRB
About Gene Symbol	



Western blot analysis of UQCRB in rat Brain lysates using Complex III Subunit 7 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using UQCRB antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.