

VDAC1/Porin Rabbit mAb


ABZ3223

Product Overview


Name	VDAC1/Porin Rabbit mAb
Accession(Primary)	P21796
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 31 kDa; Observed MW: 31 kDa
Recommended Dilution	WB-1:2000-1:20000; IHC-1:200-1:2000; ICC/IF-1:500-1:1000

Target information (P21796)


Gene Name	VDAC1
About Gene Symbol	




Western blot analysis of VDAC1 in rat kidney, mouse kidney lysates using VDAC1 antibody.



Western blot analysis of VDAC1 in 293 lysates using VDAC1 antibody.



Western blot analysis of VDAC1 in Jurkat, rat Brain, C6 lysates using VDAC1 antibody



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using VDAC1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.