

VEGFA Rabbit mAb

ABZ3224


Product Overview


Name	VEGFA Rabbit mAb
Accession(Primary)	P15692
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 27 kDa; Observed MW: 25-50 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P15692)

Gene Name VEGFA

About Gene Symbol





Western blot analysis of VEGFA in RAW264.7
lysates using VEGFA antibody.