


Villin Rabbit mAb

ABZ3227

Product Overview

Name	Villin Rabbit mAb
Accession(Primary)	P09327
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 93 kDa; Observed MW: 93 kDa
Recommended Dilution	WB-1:2000-1:5000; IHC-P-1:100-1:200

Target information (P09327)

Gene Name VIL1**About Gene Symbol**

Western blot analysis of Villin in MCF-7 lysates using Villin antibody.

