

Vinculin Rabbit mAb


ABZ3229

Product Overview


Name	Vinculin Rabbit mAb
Accession(Primary)	P18206
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 124 kDa; Observed MW: 124 kDa
Recommended Dilution	WB-1:5000-1:15000; IHC-:1:400-1:3000

Target information (P18206)

Gene Name	VCL
About Gene Symbol	VCL (Vinculin) is a structural cytoskeletal protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Vinculin in K562, C6, 3T3, HeLa lysates using Vinculin antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Vinculin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.