

XPD Rabbit mAb


ABZ3231

Product Overview


Name	XPD Rabbit mAb
Accession(Primary)	P18074
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 87 kDa; Observed MW: 80 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:200-1:500

Target information (P18074)

Gene Name	ERCC2
About Gene Symbol	



Western blot analysis of XPD in Jurkat, HeLa lysates using XPD antibody.



Immunocytochemistry analysis of XPD (green) in HeLa using XPD antibody, and DAPI (blue)