

ZNF160 Rabbit mAb

ABZ3233

Product Overview

Name	ZNF160 Rabbit mAb
Accession(Primary)	Q9HCG1
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 94 kDa; Observed MW: 74 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20

Target information (Q9HCG1)

Gene Name ZNF160**About Gene Symbol**

