

Tight Junction Protein ZO 3 Rabbit mAb

ABZ3234


Product Overview

Name	Tight Junction Protein ZO 3 Rabbit mAb
Accession(Primary)	O95049
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 101 kDa; Observed MW: 140 kDa
Recommended Dilution	WB-1:2000-1:5000

Target information (O95049)

Gene Name TJP3

About Gene Symbol



Western blot analysis of Zonula occludens protein 3 in MCF-7, Caco2 lysates using Tight Junction Protein ZO 3 antibody.

