


ATF3 Rabbit mAb

ABZ3235

Product Overview

Name	ATF3 Rabbit mAb
Accession(Primary)	P18847
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 21 kDa; Observed MW: 21 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P18847)

Gene Name ATF3**About Gene Symbol**

Western blot analysis of ATF3 in A549, MCF-7, U251, HL-60 lysates using ATF3 antibody.

