


SURF4 Rabbit mAb

ABZ3236


Product Overview

Name	SURF4 Rabbit mAb
Accession(Primary)	O15260
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 30 kDa; Observed MW: 30 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:20-1:50; FC-1:20


Target information (O15260)

Gene Name SURF4**About Gene Symbol**

Western blot analysis of SURF4 in 3T3, C2C12 lysates using Surfeit Locus Protein 4 antibody.



Western blot analysis of SURF4 in C6, Jurkat lysates using SURF4 antibody.



Immunocytochemistry analysis of SURF4 (green) in HeLa using SURF4 antibody, and DAPI (blue)