

Tissue Factor Rabbit mAb


ABZ3237

Product Overview


Name	Tissue Factor Rabbit mAb
Accession(Primary)	P13726
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 33 kDa; Observed MW: 45 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:500-1:2000; ICC/IF-1:50-1:100

Target information (P13726)


Gene Name	F3
About Gene Symbol	



Western blot analysis of Tissue Factor in mouse brain, rat brain lysates using Tissue Factor antibody.



Western blot analysis of Tissue Factor in HeLa, A549, MCF-7, U251, HL-60 lysates using Tissue Factor antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Tissue Factor antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.