


TIM 1 Rabbit mAb

ABZ3238


Product Overview

Name	TIM 1 Rabbit mAb
Accession(Primary)	Q96D42
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 39 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:2000-1:5000

Target information (Q96D42)

Gene Name HAVCR1**About Gene Symbol**

Western blot analysis of TIM 1 in Raw264.7, PC-12 lysates using Hepatitis A Virus ular Receptor 1 antibody.



Western blot analysis of TIM 1 in A549, MCF-7, U251, HL-60 lysates using TIM 1 antibody.