

Osteopontin Rabbit mAb


ABZ3239

Product Overview


Name	Osteopontin Rabbit mAb
Accession(Primary)	P10451
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 35 kDa; Observed MW: 35-100 kDa
Recommended Dilution	WB-1:500-1:4000; IHC-1:200-1:5000; FC-1:20-1:50; IHC-1:50-1:100

Target information (P10451)

Gene Name	SPP1
About Gene Symbol	



Western blot analysis of Osteopontin in HeLa, CHO-K1, C6, Jurkat lysates using Osteopontin antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Osteopontin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.