


Hex Rabbit mAb

ABZ3240

Product Overview

Name	Hex Rabbit mAb
Accession(Primary)	Q03014
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 30 kDa; Observed MW: 30 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20-1:50

Target information (Q03014)

Gene Name HHEX**About Gene Symbol**

Western blot analysis of Hex in Jurkat, HeLa lysates using Hex antibody.

