


ETS1 Rabbit mAb

ABZ3241

Product Overview

Name	ETS1 Rabbit mAb
Accession(Primary)	P14921
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 50 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:100-1:1000

Target information (P14921)

Gene Name ETS1**About Gene Symbol**

Western blot analysis of ETS1 in HeLa, A549, HL-60, U2OS lysates using ETS1 antibody.

