

PIM1 Rabbit mAb


ABZ3242

Product Overview


Name	PIM1 Rabbit mAb
Accession(Primary)	P11309
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 36 kDa; Observed MW: 36 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:20-1:50; FC-1:20

Target information (P11309)

Gene Name	PIM1
About Gene Symbol	



Western blot analysis of PIM1 in HeLa, A549, HL-60 lysates using PIM1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using PIM1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.