


IFRD1 Rabbit mAb

ABZ3243

Product Overview

Name	IFRD1 Rabbit mAb
Accession(Primary)	O00458
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 50 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200

Target information (O00458)

Gene Name IFRD1**About Gene Symbol**

Western blot analysis of IFRD1 in Ramos lysates using IFRD1 antibody.

