

Elongation Factor Ts Rabbit mAb


ABZ3244

Product Overview

Name	Elongation Factor Ts Rabbit mAb
Accession(Primary)	P43897
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 35 kDa; Observed MW: 35 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:50-1:100

Target information (P43897)

Gene Name	TSFM
About Gene Symbol	TSFM (Ts Translation Elongation Factor, Mitochondrial) is a protein kinase and signaling regulator involved in translational regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of TSFM in HeLa, A549, HL-60, U2OS lysates using Elongation Factor Ts antibody.