


SF20 Rabbit mAb

ABZ3245

Product Overview

Name	SF20 Rabbit mAb
Accession(Primary)	Q969H8
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 19 kDa; Observed MW: 19 kDa
Recommended Dilution	WB-1:500-1:1000; FC-1:20-1:50

Target information (Q969H8)

Gene Name MYDGF**About Gene Symbol**

Western blot analysis of SF20 in HeLa, Jurkat lysates using SF20 antibody.

