

FGFR3 Rabbit mAb


ABZ3246

Product Overview


Name	FGFR3 Rabbit mAb
Accession(Primary)	P22607
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 88 kDa; Observed MW: 125 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:200-1:500; FC-1:20-1:50

Target information (P22607)


Gene Name	FGFR3
About Gene Symbol	FGFR3 (Fibroblast Growth Factor Receptor 3) is a receptor tyrosine kinase involved in adaptive immunity. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of FGFR3 in HeLa lysates using FGFR3 antibody.



Immunocytochemistry analysis of FGFR3 (green) in HeLa using FGFR3 antibody, and DAPI (blue)



Immunohistochemistry analysis of paraffin-embedded Human tonsil using FGFR3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.