


RASGRP1 Rabbit mAb

ABZ3249

Product Overview

Name	RASGRP1 Rabbit mAb
Accession(Primary)	O95267
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 90 kDa; Observed MW: 90 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (O95267)

Gene Name RASGRP1**About Gene Symbol**

Western blot analysis of Rasgrp1 in A549, HL-60, U251 lysates using RASGRP1 antibody.

