


CBX4 Rabbit mAb

ABZ3250


Product Overview

Name	CBX4 Rabbit mAb
Accession(Primary)	O00257
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 61 kDa; Observed MW: 61 kDa
Recommended Dilution	WB-1:2000-1:5000

Target information (O00257)

Gene Name CBX4**About Gene Symbol**

Western blot analysis of CBX4 in C6, Raw264.7 lysates using CBX4 antibody.



Western blot analysis of CBX4 in HeLa, A549, HL-60, U251 lysates using CBX4 antibody.