

5T4 Rabbit mAb


ABZ3251

Product Overview


Name	5T4 Rabbit mAb
Accession(Primary)	Q13641
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 46 kDa; Observed MW: 78 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:50-1:100

Target information (Q13641)

Gene Name	TPBG
About Gene Symbol	TPBG (Trophoblast Glycoprotein) is a protein kinase and signaling regulator involved in Wnt/ β -catenin signaling and adaptive immunity. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of 5T4 in mouse brain, rat brain lysates using 5T4 antibody.



Western blot analysis of 5T4 in 293T, C2C12 lysates using 5T4 antibody.