

A-Raf Rabbit mAb


ABZ3252

Product Overview


Name	A-Raf Rabbit mAb
Accession(Primary)	P10398
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 68 kDa; Observed MW: 68 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:1000; ICC/IF-1:100-1:200; FC-1:50-1:100

Target information (P10398)

Gene Name	ARAF
About Gene Symbol	ARAF (A-Raf Proto-Oncogene, Serine/Threonine Kinase) is a serine/threonine protein kinase involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of A RAF in HeLa, Ramos lysates using A RAF antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using A RAF antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.