

ACADM/MCAD Rabbit mAb


ABZ3253

Product Overview


Name	ACADM/MCAD Rabbit mAb
Accession(Primary)	P11310
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:50-1:100

Target information (P11310)


Gene Name	ACADM
About Gene Symbol	ACADM (Acyl-CoA Dehydrogenase Medium Chain 2 3 5) is a metabolic dehydrogenase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mech




Immunocytochemistry analysis of ACADM (green) in 293T using ACADM antibody, and DAPI (blue).



Western blot analysis of ACADM/MCAD in rat muscle lysates using ACADM/MCAD antibody.



Western blot analysis of ACADM/MCAD in HeLa, A549, MCF-7, U251, C2C12 lysates using ACADM/MCAD antibody



Immunohistochemistry analysis of paraffin-embedded Human tonsil using ACADM/MCAD antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.