

AKAP95 Rabbit mAb


ABZ3255

Product Overview


Name	AKAP95 Rabbit mAb
Accession(Primary)	O43823
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 76 kDa; Observed MW: 95 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200

Target information (O43823)


Gene Name	AKAP8
About Gene Symbol	AKAP8 (A-Kinase Anchoring Protein 8) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Western blot analysis of AKAP 95 in PC-12, 3T3 lysates using AKAP95 antibody.



Western blot analysis of AKAP 95 in HeLa, A549, HL-60, U251 lysates using AKAP 95 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using AKAP 95 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of AKAP 95 (green) in HeLa using AKAP 95 antibody, and DAPI (blue).