

alpha Internexin Rabbit mAb


ABZ3256

Product Overview


Name	alpha Internexin Rabbit mAb
Accession(Primary)	Q16352
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 55 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:20-1:50; IF-1:20-1:50; IF-1:20; IP-1:20

Target information (Q16352)


Gene Name	INA
About Gene Symbol	




Western blot analysis of alpha Internexin in mouse brain, rat brain lysates using alpha Internexin antibody.



Western blot analysis of alpha Internexin in rat Brain lysates using alpha Internexin antibody.



Immunohistochemistry analysis of paraffin-embedded Human Brain using alpha Internexin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of alpha Internexin (green) in hela using alpha Internexin antibody, and DAPI (blue).