


Ataxin 3 Rabbit mAb

ABZ3260

Product Overview

Name	Ataxin 3 Rabbit mAb
Accession(Primary)	P54252
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 41 kDa; Observed MW: 41 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P54252)

Gene Name ATXN3**About Gene Symbol**

Western blot analysis of Ataxin 3 in HeLa, A549, HL-60, U251, U87-MG lysates using Ataxin 3 antibody.

